

IL-3 Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11858-R001

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	11858-H08B
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human IL-3 (rh IL-3; Catalog#11858-H08B; NP_000579.2; Met 1-Phe 152).
Ig Type:	Rabbit IgG
Clone ID:	001
Specificity:	Human IL-3
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.
Alternative Names:	IL-3
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000

Please Note: Optimal concentrations/dilutions should be determined by the end user.