

CD19 / B4 / CVID3 Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11880-MM02

GENERAL INFORMATION

Immunogen:	11880-H08H
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human CD19 / B4 / CVID3 (rh CD19 / B4 / CVID3; Catalog#11880-H08H; P15391-1; Met1-Lys291). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG1
Clone ID:	02
Specificity:	Human CD19 / B4 / CVID3
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.
Alternative Names:	B4

APPLICATIONS

Applications:	ELISA(Cap)
---------------	------------

RECOMMENDED CONCENTRATION

Sandwich ELISA (Capture Ab)	ELISA(Cap): 1:250-1:2000 This antibody will detect Human CD19 / B4 / CVID3 in ELISA pair set (Catalog: # SEK11880). In a sandwich ELISA, it can be used as capture antibody when paired with (Catalog: # 11880-MM15).
-----------------------------	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.