



Catalog Number: 11880-R410

EliteRmab® is a registered trademark of Sino Biological Inc.


GENERAL INFORMATION	
Immunogen:	11880-H08H
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human CD19 (rh CD19; Catalog#11880-H08H; P15391-1; Met1-Lys291).
Ig Type:	Rabbit IgG
Clone ID:	410
Specificity:	Human CD19
	No cross-reactivity in ELISA with Mouse CD19
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.
Alternative Names:	B4
APPLICATIONS	
Applications:	ELISA, IP
RECOMMENDED CONCENTRATION	
Immunoprecipitation	IP: 1-4 µL/mg of lysate
ELISA	ELISA: 1:5000-1:10000

Please Note: Optimal concentrations/dilutions should be determined by the end user.



Catalog Number: 11880-R410

EliteRmab® is a registered trademark of Sino Biological Inc.



CD19 was immunoprecipitated using:

Lane A: 0.5 mg Daudi Whole Cell Lysate

2 μ L anti-CD19 rabbit monoclonal antibody
and 15 μ L of 50 % Protein G agarose.

Primary antibody:

Anti-CD19 rabbit monoclonal antibody, at 1:100
dilution

Secondary antibody:

Dylight 800-labeled antibody to rabbit IgG
(H+L), at 1:5000 dilution

Developed using the odyssey technique.

Performed under reducing conditions.

Predicted band size: 61 kDa

Observed band size: 100 kDa