

KIT / c-KIT / CD117 Antibody (HRP), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11996-MM14-H

GENERAL INFORMATION

Immunogen:	Recombinant Human KIT / c-KIT / CD117 Protein (Catalog#11996-H08H)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human KIT / c-KIT / CD117 (rh KIT / c-KIT / CD117; Catalog#11996-H08H; P10721-2; Met1-Thr516). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography and conjugated with horseradish-peroxidase (HRP).
Ig Type:	Mouse IgG1
Clone ID:	14
Specificity:	Human KIT / c-KIT / CD117
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Protected from prolonged exposure to light. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	KIT

APPLICATIONS

Applications:	ELISA
---------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 0.1-1 µg/ml
-------	--------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.