



Catalog Number: 12042-R005


EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human GSTM2 protein (Catalog#12042-H07E)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human GSTM2 (rh GSTM2; Catalog#12042-H07E; P28161; Met1-Lys218).
Ig Type:	Rabbit IgG
Clone ID:	005
Specificity:	Human GSTM2
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	IHC-P
RECOMMENDED CONCENTRATION	
IHC-P	IHC-P: 1:100-1:500


Please Note: Optimal concentrations/dilutions should be determined by the end user.

Catalog Number: 12042-R005


EliteRmab® is a registered trademark of Sino Biological Inc.



Immunohistochemical staining of human GSTM2 in human leiomyoma with rabbit monoclonal antibody (1:200, formalin-fixed paraffin embedded sections).



Immunohistochemical staining of human GSTM2 in human heart with rabbit monoclonal antibody (1:200, formalin-fixed paraffin embedded sections).



Immunohistochemical staining of human GSTM2 in human epididymis with rabbit monoclonal antibody (1:200, formalin-fixed paraffin embedded sections). Positive staining was localized to skeletal muscle.