

CXCL3 / GRO gamma Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 12186-MM05

GENERAL INFORMATION

Immunogen:	Recombinant Human CXCL3 / GRO gamma Protein (Catalog#12186-H07E)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human CXCL3 / GRO gamma (rh CXCL3 / GRO gamma; Catalog#12186-H07E; NP_002081.2; Ala35-Asn107). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG1
Clone ID:	05
Specificity:	Human CXCL3 / GRO gamma
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human CXCL3 / GRO gamma.
--------------	---

Please Note: Optimal concentrations/dilutions should be determined by the end user.