

CD164 Antibody (FITC), Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 12260-R006-F

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	Recombinant Human CD164 Protein (Catalog#12260-H08H)
Reagents:	FITC-conjugated Rabbit monoclonal antibody
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human CD164 (rh CD164; Catalog#12260-H08H; Q04900-1; Met1-Asp162) and conjugated with FITC under optimum conditions, the unreacted FITC was removed.
Ig Type:	Rabbit IgG
Clone ID:	006
Specificity:	Human CD164
Concentration:	10 µl/Test, 0.1 mg/ml
Formulation:	PBS solution containing 0.5% BSA and 0.03%ProClin300
Storage:	This antibody can be stored at 2℃-8℃ for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze !
Alternative Names:	CD164

APPLICATIONS

Applications:	FCM
----------------------	-----


RECOMMENDED CONCENTRATION

Please Note: Optimal concentrations/dilutions should be determined by the end user.



Catalog Number: 12260-R006-F

EliteRmab® is a registered trademark of Sino Biological Inc.



Flow cytometric analysis of Human CD164 expression on human whole blood monocytes. Cells were stained with FITC-conjugated anti-Human CD164. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of viable monocytes.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.