

Carbonic Anhydrase VB / CA5B Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 12285-R008

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human Carbonic Anhydrase VB / CA5B protein (Catalog#12285-H08E)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human Carbonic Anhydrase VB / CA5B (rh Carbonic Anhydrase VB / CA5B; Catalog#12285-H08E; Q9Y2D0; Cys34-Pro317).
Ig Type:	Rabbit IgG
Clone ID:	008
Specificity:	Human Carbonic Anhydrase VB / CA5B
	Has cross-reactivity in ELISA with Human CA2 Human CA3 Human CA5A No cross-reactivity in ELISA with Human CA8 Human CA4
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	IHC-P
RECOMMENDED CONCENTRATION	
IHC-P	IHC-P: 1:50-1:200

Please Note: Optimal concentrations/dilutions should be determined by the end user.


Carbonic Anhydrase VB / CA5B Antibody, Rabbit MAb

Catalog Number: 12285-R008




Sino Biological
Biological Solution Specialist


EliteRmab® is a registered trademark of Sino Biological Inc.



Immunohistochemical staining of human CA5B in human placenta with rabbit monoclonal antibody (1:100, formalin-fixed paraffin embedded sections).



Immunohistochemical staining of human CA5B in human lung with rabbit monoclonal antibody (1:100, formalin-fixed paraffin embedded sections).



Immunohistochemical staining of human CA5B in human kidney with rabbit monoclonal antibody (1:100, formalin-fixed paraffin embedded sections).