

CA5B Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 12285-T18

GENERAL INFORMATION

Immunogen:	Recombinant Human CA5B Protein (Catalog#12285-H08E)
Preparation	Produced in rabbits immunized with purified, recombinant Human CA5B (rh CA5B; Catalog#12285-H08E; Q9Y2D0; Cys34-Pro317). CA5B specific IgG was purified by Human CA5B affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human CA5B
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	CA-VB,MGC39962

APPLICATIONS

Applications:	ELISA, ICC/IF
---------------	---------------

RECOMMENDED CONCENTRATION

ICC/IF	ICC/IF: 1:300-1:10000
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human CA5B.


Please Note: Optimal concentrations/dilutions should be determined by the end user.

CA5B Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 12285-T18



Immunofluorescence staining of CA5B in MCF7 cells. Cells were fixed with 4% PFA, permeabilized with 0.1% Triton X-100 in PBS, blocked with 10% serum, and incubated with rabbit anti-Human CA5B polyclonal antibody (dilution ratio 1:1000) at 4°C overnight. Then cells were stained with the Alexa Fluor®488-conjugated Goat Anti-rabbit IgG secondary antibody (green). Positive staining was localized to Cytoplasm.