

AK4 / Adenylate Kinase 4 / AK3L1 Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 12406-RP01

GENERAL INFORMATION	
Immunogen:	Recombinant Human AK4 / Adenylate Kinase 4 / AK3L1 protein (Catalog#12406-H20B)
Preparation	Produced in rabbits immunized with purified, recombinant Human AK4 / Adenylate Kinase 4 / AK3L1 (rh AK4; Catalog#12406-H20B; P27144-1; Ala2-Tyr223). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human AK4 / Adenylate Kinase 4 / AK3L1
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human AK4 / Adenylate Kinase 4 / AK3L1.

Please Note: Optimal concentrations/dilutions should be determined by the end user.