

DDOST / OST48 Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 12463-MM02

GENERAL INFORMATION

Immunogen:	Recombinant Human DDOST / OST48 protein (Catalog#12463-HNCE)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human DDOST / OST48 (rh DDOST / OST48; Catalog#12463-HNCE; P39656-1; Ser43-Pro427). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG2a
Clone ID:	02
Specificity:	Human DDOST / OST48
	No cross-reactivity in ELISA with E.coli cell lysate
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications: ELISA

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human DDOST.
-------	---

Please Note: Optimal concentrations/dilutions should be determined by the end user.