

PHPT1 Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 12473-MM11

GENERAL INFORMATION

Immunogen:	Recombinant Human PHPT1 Protein (Catalog#12473-H07E)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human PHPT1 (rh PHPT1; Catalog#12473-H07E; Q9NRX4-1; Ala2-Tyr125). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG1
Clone ID:	11
Specificity:	Human PHPT1
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA(Cap)
----------------------	------------

RECOMMENDED CONCENTRATION

Sandwich ELISA (Capture Ab)	ELISA(Cap): 1:250-1:2000 This antibody will detect Human PHPT1 in ELISA pair set (Catalog: # SEK12473). In a sandwich ELISA, it can be used as capture antibody when paired with (Catalog: # 12473-MM01).
--	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.