

FABP6 / ILBP / I-BABP Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 12579-RP02

GENERAL INFORMATION	
Immunogen:	Recombinant Human FABP6 / ILBP / I-BABP protein (Catalog#12579-HNAE)
Preparation	Produced in rabbits immunized with purified, recombinant Human FABP6 / ILBP / I-BABP (rh FABP6 / ILBP / I-BABP; Catalog#12579-HNAE; AAH22489.1; Met 1-Ala 128). FABP6 / ILBP / I-BABP specific IgG was purified by Human FABP6 / ILBP / I-BABP affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human FABP6 / ILBP / I-BABP
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human FABP6.

Please Note: Optimal concentrations/dilutions should be determined by the end user.