

IGF-2 / IGF-II Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 13032-T16

GENERAL INFORMATION

Immunogen:	Recombinant Human IGF-2 / IGF-II Protein (Catalog#13032-HNAY)
Preparation	Produced in rabbits immunized with purified, recombinant Human IGF-2 / IGF-II (rh IGF-2 / IGF-II; Catalog#13032-HNAY; P01344-1; Ala25-Glu91). IGF-2 / IGF-II specific IgG was purified by Human IGF-2 / IGF-II affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human IGF-2 / IGF-II
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	IGF2

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human IGF-2 / IGF-II.
--------------	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.