

IL22 / Interleukin 22 Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 13059-MM24

GENERAL INFORMATION

Immunogen:	Recombinant Human IL22 / Interleukin 22 Protein (Catalog#13059-HNAE)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human IL22 / Interleukin 22 (rh IL22 / Interleukin 22; Catalog#13059-HNAE; NP_065386.1; Ala34-Ile179). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG1
Clone ID:	24
Specificity:	Human IL22 / Interleukin 22
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA(Cap)
----------------------	------------

RECOMMENDED CONCENTRATION

Sandwich ELISA (Capture Ab)	ELISA(Cap): 1:250-1:2000 This antibody will detect Human IL22 / Interleukin 22 in ELISA pair set (Catalog: # SEK13059). In a sandwich ELISA, it can be used as capture antibody when paired with (Catalog: # 13059-R034).
------------------------------------	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.