

# Coagulation Factor III / Tissue Factor / CD142 Antibody (APC), Mouse MAb



Sino Biological  
Biological Solution Specialist

Catalog Number: 13133-MM05-A

## GENERAL INFORMATION

|                           |                                                                                                                                                                                                                                                                                                                                                               |
|---------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Immunogen:</b>         | Recombinant Human F3 / Tissue factor / CD142 protein (Catalog#13133-H08H)                                                                                                                                                                                                                                                                                     |
| <b>Reagents:</b>          | APC-conjugated Mouse monoclonal antibody                                                                                                                                                                                                                                                                                                                      |
| <b>Preparation</b>        | This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human F3 / Tissue factor / CD142 (rh F3 / Tissue factor / CD142; Catalog#13133-H08H; P13726; Met 1-Glu 251) and conjugated with APC under optimum conditions, the unreacted APC was removed. |
| <b>Ig Type:</b>           | Mouse IgG2a                                                                                                                                                                                                                                                                                                                                                   |
| <b>Clone ID:</b>          | 05                                                                                                                                                                                                                                                                                                                                                            |
| <b>Specificity:</b>       | Human F3 / Tissue factor / CD142                                                                                                                                                                                                                                                                                                                              |
| <b>Concentration:</b>     | 1 µl/Test, 0.1 mg/ml                                                                                                                                                                                                                                                                                                                                          |
| <b>Formulation:</b>       | PBS solution containing 0.5% BSA and 0.03%ProClin300                                                                                                                                                                                                                                                                                                          |
| <b>Storage:</b>           | This antibody can be stored at 2°C-8°C for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze !                                                                                                                                                                                                     |
| <b>Alternative Names:</b> | F3                                                                                                                                                                                                                                                                                                                                                            |

## APPLICATIONS

|                      |     |
|----------------------|-----|
| <b>Applications:</b> | FCM |
|----------------------|-----|

## RECOMMENDED CONCENTRATION


*Please Note: Optimal concentrations/dilutions should be determined by the end user.*

# Coagulation Factor III / Tissue Factor / CD142 Antibody (APC), Mouse MAb



Sino Biological  
Biological Solution Specialist

Catalog Number: 13133-MM05-A



Flow cytometric analysis of F3 (CD142) expression on human A431 cells.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system.  
Please refer to [www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html](http://www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html) for technical protocols.

**For Research Use Only. Not for use in diagnostic or therapeutic procedures.**

**Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)**

**Website: <http://www.sinobiological.com>**