

Protein HC/Bikunin Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 13141-MM14

GENERAL INFORMATION

Immunogen:	Recombinant Human Protein HC/Bikunin Protein (Catalog#13141-H08H1)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human Protein HC/Bikunin (rh Protein HC/Bikunin; Catalog#13141-H08H1; NP_001624.1; Met1-Val203). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG1
Clone ID:	14
Specificity:	Human Protein HC/Bikunin
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
---------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human Protein HC/Bikunin.
-------	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.