

CD3 zeta / CD247 Antibody (PE), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 13142-MM06-P

GENERAL INFORMATION

Immunogen:	Recombinant Human CD3 zeta / CD247 Protein (Catalog#13142-H09E)
Reagents:	PE-conjugated Mouse MAb antibody
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human CD3 zeta / CD247 (rh CD3 zeta / CD247; Catalog#13142-H09E; P20963-1; Arg52-Arg164) and conjugated with PE under optimum conditions, the unreacted PE was removed.
Ig Type:	Mouse IgG2a
Clone ID:	06
Specificity:	Human CD3 zeta / CD247
Concentration:	5 µl/Test, 0.1 mg/ml
Formulation:	PBS solution containing 0.5% BSA and 0.03%ProClin300
Storage:	This antibody can be stored at 2°C-8°C for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze !

APPLICATIONS

Applications:	FCM
----------------------	-----

RECOMMENDED CONCENTRATION


Please Note: Optimal concentrations/dilutions should be determined by the end user.

CD3 zeta / CD247 Antibody (PE), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 13142-MM06-P



Flow cytometric analysis of Human CD3 Zeta expression on human whole blood lymphocytes. The cells were treated according to manufacturer's manual (BD Pharmingen™ Cat. No. 554714), stained with PE-conjugated anti-Human CD3 Zeta and FITC-conjugated anti-Human CD3. The dot plots were derived from gated events with the forward and side light-scatter characteristics of viable lymphocytes.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>