

CCDC134 Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 13168-MM04

GENERAL INFORMATION	
Immunogen:	Recombinant Human CCDC134 protein (Catalog#13168-H08H)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human CCDC134 (rh CCDC134; Catalog#13168-H08H; Q9H6E4; Met 1-Leu 229). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG1
Clone ID:	04
Specificity:	Human CCDC134
	No cross-reactivity in ELISA with Human cell lysate (293 cell line)
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human CCDC134.

Please Note: Optimal concentrations/dilutions should be determined by the end user.