

MAG Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 13186-T16

GENERAL INFORMATION

Immunogen:	Recombinant Human MAG Protein (Catalog#13186-H08H)
Preparation	Produced in rabbits immunized with purified, recombinant Human MAG (rh MAG; Catalog#13186-H08H; P20916; Met1-Pro516). MAG specific IgG was purified by Human MAG affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human MAG
Formulation:	PBS, pH7.0 with 0.03% Proclin300
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.
Alternative Names:	MAG

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human MAG.
--------------	---

Please Note: Optimal concentrations/dilutions should be determined by the end user.