

APLP1 / Amyloid-like protein 1 Antibody, Rabbit PAb, Antigen Affinity Purified

Catalog Number: 13716-T16



Sino Biological
Biological Solution Specialist

GENERAL INFORMATION

Immunogen:	Recombinant Human APLP1 / Amyloid-like protein 1 protein (Catalog#13716-H08H)
Preparation	Produced in rabbits immunized with purified, recombinant Human APLP1 / Amyloid-like protein 1 (rh APLP1 / Amyloid-like protein 1; Catalog#13716-H08H; P51693-1; Met1-Glu580). APLP1 / Amyloid-like protein 1 specific IgG was purified by Human APLP1 / Amyloid-like protein 1 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human APLP1 / Amyloid-like protein 1
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human APLP1 / Amyloid-like protein 1.
--------------	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.