

AMIGO2 Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 13735-RP01

GENERAL INFORMATION

Immunogen:	Recombinant Human AMIGO2 protein (Catalog#13735-H08H)
Preparation	Produced in rabbits immunized with purified, recombinant Human AMIGO2 (rh AMIGO2; Catalog#13735-H08H; NP_862830.1; Met1-Thr398). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human AMIGO2
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
---------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human AMIGO2.
-------	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.