

BPI Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 13907-RP01

GENERAL INFORMATION	
Immunogen:	Recombinant Human BPI protein (Catalog#13907-H08H)
Preparation	Produced in rabbits immunized with purified, recombinant Human BPI (rh BPI; Catalog#13907-H08H; P17213; Met1-Lys487). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human BPI
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human BPI.

Please Note: Optimal concentrations/dilutions should be determined by the end user.