

Coagulation Factor VIII / FVIII / F8 (Light Chain) Antibody (Biotin), Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 13909-R402-B

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	Recombinant Human Coagulation Factor VIII / FVIII / F8 protein
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human Coagulation Factor VIII / FVIII / F8 protein and then biotinylated.
Ig Type:	Rabbit IgG
Clone ID:	402
Specificity:	Human Coagulation Factor VIII / FVIII / F8 Light Chain
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 0.1-1 µg/ml
--------------	--------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.