

BPIFA2 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 13915-T16

GENERAL INFORMATION

Immunogen:	Recombinant Human BPIFA2 Protein (Catalog#13915-H08H)
Preparation	Produced in rabbits immunized with purified, recombinant Human BPIFA2 (rh BPIFA2; Catalog#13915-H08H; Q96DR5; Met1-Ile249). BPIFA2 specific IgG was purified by Human BPIFA2 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human BPIFA2
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
---------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000
-------	-----------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.
Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)
Website: <http://www.sinobiological.com>