

LSP1 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 14178-T56

GENERAL INFORMATION

Immunogen:	Recombinant Human LSP1 Protein (Catalog#14178-H08H)
Preparation	Produced in rabbits immunized with purified, recombinant Human LSP1 (rh LSP1; Catalog#14178-H08H; AAH01785.1; Met1-Pro339). LSP1 specific IgG was purified by Human LSP1 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human LSP1
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	WB,ELISA,IHC-P
----------------------	----------------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:1000-1:4000
Western Blot	WB: 1:500-1:2000
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human LSP1.


Please Note: Optimal concentrations/dilutions should be determined by the end user.

LSP1 Antibody, Rabbit PAb, Antigen Affinity Purified




Sino Biological
Biological Solution Specialist


Catalog Number: 14178-T56



Immunohistochemical staining of human LSP1 in human tonsil with rabbit polyclonal antibody (1:2000, formalin-fixed paraffin embedded sections).



Immunohistochemical staining of human LSP1 in human spleen with rabbit polyclonal antibody (1:2000, formalin-fixed paraffin embedded sections).



Anti-S3-LSP1 rabbit polyclonal antibody at 1:500 dilution

Lane A: Daudi Whole Cell Lysate

Lysates/proteins at 30 µg per lane.

Secondary

Goat Anti-Rabbit IgG H&L (Dylight800) at 1/10000 dilution.

Developed using the Odyssey technique.
Performed under reducing conditions.

Predicted band size: 37 kDa

Observed band size: 50 kDa