

CD177 / PRV-1 / PRV1 Antibody (PE), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 14501-MM05-P

GENERAL INFORMATION

Immunogen:	Recombinant Human CD177 / PRV-1 / PRV1 Protein (Catalog#14501-H08H)
Reagents:	PE-conjugated Mouse monoclonal antibody
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human CD177 / PRV-1 / PRV1 (rh CD177 / PRV-1 / PRV1; Catalog#14501-H08H; AAH29167.1; Met1-Gly407) and conjugated with PE under optimum conditions, the unreacted PE was removed.
Ig Type:	Mouse IgG1
Clone ID:	05
Specificity:	Human CD177 / PRV-1 / PRV1
Concentration:	5 µl/Test, 0.1 mg/ml
Formulation:	PBS solution containing 0.5% BSA and 0.09% sodium azide
Storage:	This antibody can be stored at 2°C-8°C for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
Alternative Names:	CD177

APPLICATIONS

Applications:	FCM
----------------------	-----

RECOMMENDED CONCENTRATION


Please Note: Optimal concentrations/dilutions should be determined by the end user.

CD177 / PRV-1 / PRV1 Antibody (PE), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 14501-MM05-P



Flow cytometric analysis of Human CD177 expression on human whole blood granulocytes. Cells were stained with PE-conjugated anti-Human CD177. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of viable granulocytes.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.