

EIF3K Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 14553-MM01

GENERAL INFORMATION	
Immunogen:	Recombinant Human EIF3K Protein (Catalog#14553-H20B)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human EIF3K (rh EIF3K; Catalog#14553-H20B; Q9UBQ5; Ala2-Gln218). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG1
Clone ID:	01
Specificity:	Human EIF3K
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	ARG134,EIF3-p28,EIF3S12,HSPC029,M9,MSTP001,PLAC-24,PLAC24,PRO1474,PTD001
APPLICATIONS	
Applications:	ICC/IF
RECOMMENDED CONCENTRATION	
ICC/IF	ICC/IF: 1:20-1:100


Please Note: Optimal concentrations/dilutions should be determined by the end user.

EIF3K Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 14553-MM01



Immunofluorescence staining of Human EIF3K in HeLa cells. Cells were fixed with 4% PFA, permeabilized with 0.3% Triton X-100 in PBS, blocked with 10% serum, and incubated with mouse anti-Human EIF3K monoclonal antibody (1:60) at 4°C overnight. Then cells were stained with the Alexa Fluor® 488-conjugated Goat Anti-mouse IgG secondary antibody (green). Positive staining was localized to cytoplasm and nucleus.