

EIF5A2 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 14927-T48

GENERAL INFORMATION

Immunogen:	Recombinant Human EIF5A2 Protein (Catalog#14927-H07E)
Preparation	Produced in rabbits immunized with purified, recombinant Human EIF5A2 (rh EIF5A2; Catalog#14927-H07E; Q9GZV4; Met1-Lys153). EIF5A2 specific IgG was purified by Human EIF5A2 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human EIF5A2
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	EIF-5A2, eIF5AII

APPLICATIONS

Applications:	WB, ELISA
----------------------	-----------

RECOMMENDED CONCENTRATION

Western Blot	WB: 1:500-1:2000
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human EIF5A2.


Please Note: Optimal concentrations/dilutions should be determined by the end user.

EIF5A2 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 14927-T48



Anti-EIF5A2 rabbit polyclonal antibody at
1:500 dilution

Lane A: HeLa Whole Cell Lysate

Lane B: A431 Whole Cell Lysate

Lane C: Caco-2 Whole Cell Lysate

Lane D: HepG2 Whole Cell Lysate

Lysates/proteins at 30 µg per lane.

Secondary

Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000
dilution.

Developed using the ECL technique.

Performed under reducing conditions.

Predicted band size:14 kDa

Observed band size:18 kDa

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>