

ADPRH Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 15393-T52

GENERAL INFORMATION	
Immunogen:	Recombinant Human ADPRH Protein (Catalog#15393-H07E)
Preparation	Produced in rabbits immunized with purified, recombinant Human ADPRH (rh ADPRH; Catalog#15393-H07E; NP_001116.1; Met1-Leu351). ADPRH specific IgG was purified by Human ADPRH affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human ADPRH
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	WB, ELISA, IP
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000
Western Blot	WB: 1:500-1:2000
IP	IP:1-5 µL/mg of lysate


Please Note: Optimal concentrations/dilutions should be determined by the end user.

ADPRH Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 15393-T52



Anti-ADPRH rabbit polyclonal antibody at 1:500 dilution
Lane A: HeLa whole cell lysate
Lane B: K562 whole cell lysate

Lysates/proteins at 30 µg per lane.

Secondary


Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution.

Developed using the ECL technique.

Performed under reducing conditions.

Predicted band size: 40 kDa

Observed band size: 39 kDa



ADPRH was immunoprecipitated using:
Lane A: 0.5 mg K562 Whole Cell Lysate

4 µL anti-ADPRH rabbit polyclonal antibody and 60 µg of
Immunomagnetic beads Protein A/G.

Primary antibody:

Anti-ADPRH rabbit polyclonal antibody, at 1:100 dilution

Secondary antibody:

Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution

Developed using the ECL technique.

Performed under reducing conditions.

Predicted band size: 38 kDa

Observed band size: 38 kDa