

IL-5/IL5 / Interleukin 5 Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 15673-R001

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human IL-5/IL5 / Interleukin 5 Protein (Catalog#15673-HNCE)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human IL-5/IL5 / Interleukin 5 (rh IL-5/IL5 / Interleukin 5; Catalog#15673-HNCE; P05113; Ile20-Ser134).
Ig Type:	Rabbit IgG
Clone ID:	001
Specificity:	Human IL-5/IL5 / Interleukin 5
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA(Cap)
RECOMMENDED CONCENTRATION	
Sandwich ELISA (Capture Ab)	ELISA(Cap): 1:250-1:2000 This antibody will detect Human IL-5/IL5 / Interleukin 5 in ELISA pair set (Catalog: # SEKA15673). In a sandwich ELISA, it can be used as capture antibody when paired with (Catalog: # 15673-R013).

Please Note: Optimal concentrations/dilutions should be determined by the end user.