

Anti-Human IgG1 F(ab)2 Antibody (HRP), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 16008-MM02-H

GENERAL INFORMATION

Preparation:	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified. The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography and conjugated with horseradish-peroxidase (HRP).
Ig Type:	Mouse IgG1
Clone ID:	02
Specificity:	Recognize the Human IgG1 full-length antibody and F(ab)2 in ELISA, don't recognize the Human IgG2, Human IgG4, and Human IgG Fc region.
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Protected from prolonged exposure to light. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 0.1-1 µg/mL
--------------	--------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.