

LY75 / DEC-205 / CD205 Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 16490-MM16

GENERAL INFORMATION	
Immunogen:	Recombinant Human LY75 / DEC-205 / CD205 Protein (Catalog#16490-H08H)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human LY75 / DEC-205 / CD205 (rh LY75 / DEC-205 / CD205; Catalog#16490-H08H; NP_002340.2;Ala32-Gly1664). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG1
Clone ID:	16
Specificity:	Human LY75 / DEC-205 / CD205
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	LY75
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human LY75 / DEC-205 / CD205.

Please Note: Optimal concentrations/dilutions should be determined by the end user.