

AMH Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 17500-MM14

GENERAL INFORMATION	
Immunogen:	Recombinant Human AMH Protein (Catalog#17500-H07H)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human AMH (rh AMH; Catalog#17500-H07H; NP_000470.2; Leu25-Arg560). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG1
Clone ID:	14
Specificity:	Human AMH
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	MIF,MIS
APPLICATIONS	
Applications:	IHC-P
RECOMMENDED CONCENTRATION	
IHC-P	IHC-P: 1:20-1:50


Please Note: Optimal concentrations/dilutions should be determined by the end user.

AMH Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 17500-MM14



Immunohistochemical staining of human AMH in human testis with mouse monoclonal antibody at 1:30 dilution, formalin-fixed paraffin embedded sections.