

Myosin light chain 3 / MYL3 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 200249-T08

GENERAL INFORMATION

Immunogen:	E. coli-derived Human Myosin light chain 3 / MYL3 fragment
Preparation	Produced in rabbits immunized with E. coli-derived Human Myosin light chain 3 / MYL3 fragment, and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human Myosin light chain 3 / MYL3
Formulation:	PBS, pH7.0 with 0.03% Proclin300
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.
Alternative Names:	MYL3

APPLICATIONS

Applications:	IHC-P
---------------	-------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:200-1:1000
-------	---------------------


Please Note: Optimal concentrations/dilutions should be determined by the end user.

Myosin light chain 3 / MYL3 Antibody, Rabbit PAb, Antigen Affinity Purified




Sino Biological
Biological Solution Specialist

Catalog Number: 200249-T08



Immunohistochemical staining of human MYL3 in human skeletal muscle with rabbit polyclonal antibody at 1:500 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of human MYL3 in cynomolgus heart with rabbit polyclonal antibody at 1:500 dilution, formalin-fixed paraffin embedded sections.