

# Erbin/ERBB2IP Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological  
Biological Solution Specialist

Catalog Number: 200335-T40

## GENERAL INFORMATION

|                    |                                                                                                                                                                                                                                 |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Immunogen:         | E. coli-derived Human Erbin/ERBB2IP fragment                                                                                                                                                                                    |
| Preparation        | Produced in rabbits immunized with E. coli-derived Human Erbin/ERBB2IP fragment, and purified by antigen affinity chromatography.                                                                                               |
| Ig Type:           | Rabbit IgG                                                                                                                                                                                                                      |
| Specificity:       | Human Erbin/ERBB2IP                                                                                                                                                                                                             |
| Formulation:       | PBS, pH7.0 with 0.03% Proclin300                                                                                                                                                                                                |
| Storage:           | This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles. |
| Alternative Names: | ERBIN                                                                                                                                                                                                                           |

## APPLICATIONS

|               |           |
|---------------|-----------|
| Applications: | WB, IHC-P |
|---------------|-----------|

## RECOMMENDED CONCENTRATION

|              |                     |
|--------------|---------------------|
| IHC-P        | IHC-P: 1:500-1:2000 |
| Western Blot | WB: 1:500-1:2000    |


***Please Note: Optimal concentrations/dilutions should be determined by the end user.***

# Erbin/ERBB2IP Antibody, Rabbit PAb, Antigen Affinity Purified




Sino Biological  
Biological Solution Specialist

Catalog Number: 200335-T40



Immunohistochemical staining of human ERBB2IP in human small intestine with rabbit polyclonal antibody at 1:1000 dilution, formalin-fixed paraffin embedded sections.



Anti-ERBB2IP rabbit polyclonal antibody at 1:500 dilution

Lane A: Jurkat Whole Cell Lysate

Lysates/proteins at 30 µg per lane.

Secondary

Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution.

Developed using the ECL technique.  
Performed under reducing conditions.

Predicted band size: 158 kDa