

HEXB Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 200592-T40

GENERAL INFORMATION

Immunogen:	E. coli-derived Human HEXB fragment
Preparation	Produced in rabbits immunized with E. coli-derived Human HEXB fragment, and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human HEXB
Formulation:	PBS, pH7.0 with 0.03% Proclin300
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.
Alternative Names:	HEXB

APPLICATIONS

Applications:	WB, IHC-P
---------------	-----------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:50-1:200
Western Blot	WB: 1:500-1:2000


Please Note: Optimal concentrations/dilutions should be determined by the end user.

HEXB Antibody, Rabbit PAb, Antigen Affinity Purified


Catalog Number: 200592-T40



Sino Biological
Biological Solution Specialist



Immunohistochemical staining of human HEXB in human epididymis with rabbit polyclonal antibody at 1:100 dilution, formalin-fixed paraffin embedded sections.



Anti-HEXB rabbit polyclonal antibody at 1:500 dilution

Lane A: Jurkat Whole Cell Lysate

Lane B: HeLa Whole Cell Lysate

Lane C: HepG2 Whole Cell Lysate

Lysates/proteins at 30 µg per lane.

Secondary

Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution.

Developed using the ECL technique.

Performed under reducing conditions.

Predicted band size: 63 kDa

Observed band size: 63 kDa

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>