

HIF-1 alpha/HIF1A Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 200595-T02

GENERAL INFORMATION

Immunogen:	E. coli-derived Human HIF-1 alpha/HIF1A fragment
Preparation	Produced in rabbits immunized with E. coli-derived Human HIF-1 alpha/HIF1A fragment, and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human HIF-1 alpha/HIF1A
Formulation:	PBS, pH7.0 with 0.03% Proclin300
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	bHLHe78,HIF-1A,HIF-1alpha,HIF1,HIF1-ALPHA,MOP1,PASD8

APPLICATIONS

Applications:	ICC/IF
---------------	--------

RECOMMENDED CONCENTRATION

ICC/IF	ICC/IF: 1:100-1:500
--------	---------------------


Please Note: Optimal concentrations/dilutions should be determined by the end user.

HIF-1 alpha/HIF1A Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 200595-T02



Immunofluorescence staining of HIF1A in HeLa(500uM CoCl₂ for 24hr) cells. Cells were fixed with 4% PFA, permeabilized with 0.1% Triton X-100 in PBS, blocked with 10% serum, and incubated with rabbit anti-Human HIF1A polyclonal antibody (dilution ratio 1:200) at 4℃ overnight. Then cells were stained with the Alexa Fluor®488-conjugated Goat Anti-rabbit IgG secondary antibody (green). Positive staining was localized to Nucleus.