

BMF Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 200740-T10

GENERAL INFORMATION	
Immunogen:	E. coli-derived Human BMF fragment
Preparation	Produced in rabbits immunized with E. coli-derived Human BMF fragment, and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human BMF
Formulation:	PBS, pH7.0 with 0.03% Proclin300
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	BMF
APPLICATIONS	
Applications:	IHC-P, ICC/IF
RECOMMENDED CONCENTRATION	
IHC-P	IHC-P: 1:50-1:200
ICC/IF	ICC/IF: 1:100-1:500


Please Note: Optimal concentrations/dilutions should be determined by the end user.

BMF Antibody, Rabbit PAb, Antigen Affinity Purified


Catalog Number: 200740-T10



Sino Biological
Biological Solution Specialist



Immunohistochemical staining of human BMF in human brain with rabbit polyclonal antibody at 1:100 dilution, formalin-fixed paraffin embedded sections.



Immunofluorescence staining of BMF in HeLa cells. Cells were fixed with 4% PFA, permeabilized with 0.1% Triton X-100 in PBS, blocked with 10% serum, and incubated with rabbit anti-Human BMF polyclonal antibody (dilution ratio 1:200) at 4°C overnight. Then cells were stained with the Alexa Fluor®488-conjugated Goat Anti-rabbit IgG secondary antibody (green). Positive staining was localized to Nucleus and Cytoplasm.