

**intestinal alkaline phosphatase/ALPI
Antibody, Rabbit PAb, Antigen
Affinity Purified**

Catalog Number: 203001-T08



Sino Biological
Biological Solution Specialist

GENERAL INFORMATION

Immunogen:	E. coli-derived Human intestinal alkaline phosphatase/ALPI fragment
Preparation	Produced in rabbits immunized with E. coli-derived Human intestinal alkaline phosphatase/ALPI fragment, and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human intestinal alkaline phosphatase/ALPI
Formulation:	PBS, pH7.0 with 0.03% Proclin300
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	IHC-P
----------------------	-------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:50-1:200
--------------	-------------------


Please Note: Optimal concentrations/dilutions should be determined by the end user.

intestinal alkaline phosphatase/ALPI
Antibody, Rabbit PAb, Antigen
Affinity Purified

Catalog Number: 203001-T08



Sino Biological
Biological Solution Specialist



Immunohistochemical staining of human ALPI in human placenta with rabbit polyclonal antibody at 1:100 dilution, formalin-fixed paraffin embedded sections.