

ARID3A Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 203510-T36

GENERAL INFORMATION

Immunogen:	E. coli-derived Human ARID3A fragment
Preparation	Produced in rabbits immunized with E. coli-derived Human ARID3A fragment, and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human ARID3A
Formulation:	PBS, pH7.0 with 0.03% Proclin300
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	WB, IP
---------------	--------

RECOMMENDED CONCENTRATION

Western Blot	WB: 1:1000-1:5000
Immunoprecipitation	IP: 0.5-1µL/mg of lysate


Please Note: Optimal concentrations/dilutions should be determined by the end user.

ARID3A Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 203510-T36



Anti-ARID3A rabbit polyclonal antibody at 1:2000 dilution

Lane A: 293T Whole Cell Lysate
Lane B: HepG2 Whole Cell Lysate
Lane C: K-562 Whole Cell Lysate


Lysates/proteins at 30 µg per lane.

Secondary

Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution.

Developed using the ECL technique.
Performed under reducing conditions.

Predicted band size: 63 kDa



ARID3A was immunoprecipitated using:
Lane A: 0.5 mg HepG2 Whole Cell Lysate

4 µL anti-ARID3A rabbit polyclonal antibody
and 60 µg of Immunomagnetic beads Protein A/G.

Primary antibody:

Anti-ARID3A rabbit polyclonal antibody, at 1:100 dilution

Secondary antibody:

Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution

Developed using the ECL technique.
Performed under reducing conditions.

Predicted band size: 63 kDa
Observed band size: 80 kDa