

Annexin A1/ANXA1 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 203969-T44

GENERAL INFORMATION

Immunogen:	E. coli-derived Human Annexin A1/ANXA1 fragment
Preparation	Produced in rabbits immunized with E. coli-derived Human Annexin A1/ANXA1 fragment, and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human Annexin A1/ANXA1
Formulation:	PBS, pH7.0 with 0.03% Proclin300
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	WB, IHC-P, IP
----------------------	---------------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:5000-1:20000
Western Blot	WB: 1:500-1:2000
Immunoprecipitation	IP: 1-5 µL/mg of lysate


Please Note: Optimal concentrations/dilutions should be determined by the end user.

Annexin A1/ANXA1 Antibody, Rabbit PAb, Antigen Affinity Purified




Sino Biological
Biological Solution Specialist


Catalog Number: 203969-T44



Immunohistochemical staining of human ANXA1 in human esophagus with rabbit polyclonal antibody at 1:10000 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of human ANXA1 in human prostate with rabbit polyclonal antibody at 1:10000 dilution, formalin-fixed paraffin embedded sections.



Anti-ANXA1 rabbit polyclonal antibody at 1:500 dilution

Lane A: A431 Whole Cell Lysate

Lane B: H1299 Whole Cell Lysate

Lysates/proteins at 30 µg per lane.


Secondary

Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution.

Developed using the ECL technique.

Performed under reducing conditions.

Predicted band size:39 kDa



ANXA1 was immunoprecipitated using:

Lane A:0.5 mg A431 Whole Cell Lysate

Lane B:0.5 mg H1299 Whole Cell Lysate

4 µL anti-ANXA1 rabbit polyclonal antibody and 60 µg of Immunomagnetic beads Protein A/G.

Primary antibody:

Anti-ANXA1 rabbit polyclonal antibody,at 1:100 dilution

Secondary antibody:

Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution

Developed using the ECL technique.

Performed under reducing conditions.

Predicted band size: 39 kDa

Observed band size :37 kDa

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>