

FTSJ2 / MRM2 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 205440-T32

GENERAL INFORMATION

Immunogen:	E. coli-derived Human FTSJ2 / MRM2 fragment
Preparation	Produced in rabbits immunized with E. coli-derived Human FTSJ2 / MRM2 fragment, and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human FTSJ2 / MRM2
Formulation:	PBS, pH7.0 with 0.03% Proclin300
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.
Alternative Names:	MRM2

APPLICATIONS

Applications:	WB
---------------	----

RECOMMENDED CONCENTRATION

Western Blot	WB: 1:500-1:2000
--------------	------------------


Please Note: Optimal concentrations/dilutions should be determined by the end user.

FTSJ2 / MRM2 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 205440-T32



Anti-FTSJ2 rabbit polyclonal antibody at 1:500 dilution Lane A: U2OS Whole Cell Lysate Lane B: A431 Whole Cell Lysate Lane C: K562 Whole Cell Lysate Lane D: U-251 MG Whole Cell Lysate Lysates/proteins at 30 µg per lane. Secondary Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution. Developed using the ECL technique. Performed under reducing conditions. Predicted band size: 27 kDa Observed band size: 34 kDa