

CDC42EP2 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 207569-T08

GENERAL INFORMATION

Immunogen:	E. coli-derived Human CDC42EP2 fragment
Preparation	Produced in rabbits immunized with E. coli-derived Human CDC42EP2 fragment, and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human CDC42EP2
Formulation:	PBS, pH7.0 with 0.03% Proclin300
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.
Alternative Names:	CDC42EP2

APPLICATIONS

Applications:	IHC-P
---------------	-------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:500-1:2000
-------	---------------------


Please Note: Optimal concentrations/dilutions should be determined by the end user.

CDC42EP2 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 207569-T08



Immunohistochemical staining of human
CDC42EP2 in human kidney with rabbit
polyclonal antibody at 1:1000 dilution,
formalin-fixed paraffin embedded sections.