

# MAST4 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological  
Biological Solution Specialist

Catalog Number: 207682-T08

## GENERAL INFORMATION

|                    |                                                                                                                                                                                                                                                    |
|--------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Immunogen:         | E. coli-derived Human MAST4 fragment                                                                                                                                                                                                               |
| Preparation        | Produced in rabbits immunized with E. coli-derived Human MAST4 fragment, and purified by antigen affinity chromatography.                                                                                                                          |
| Ig Type:           | Rabbit IgG                                                                                                                                                                                                                                         |
| Specificity:       | Human MAST4                                                                                                                                                                                                                                        |
| Formulation:       | 0.2 µm filtered solution in PBS                                                                                                                                                                                                                    |
| Storage:           | This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles. |
| Alternative Names: | MAST4                                                                                                                                                                                                                                              |

## APPLICATIONS

|               |       |
|---------------|-------|
| Applications: | IHC-P |
|---------------|-------|

## RECOMMENDED CONCENTRATION

|       |                   |
|-------|-------------------|
| IHC-P | IHC-P: 1:50-1:200 |
|-------|-------------------|


**Please Note: Optimal concentrations/dilutions should be determined by the end user.**

# MAST4 Antibody, Rabbit PAb, Antigen Affinity Purified

Catalog Number: 207682-T08



Sino Biological  
Biological Solution Specialist



Immunohistochemical staining of human MAST4 in human urinary bladder with rabbit polyclonal antibody at 1:100 dilution, formalin-fixed paraffin embedded sections.