

CALML5 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 207968-T08

GENERAL INFORMATION

Immunogen:	E. coli-derived Human CALML5 fragment
Preparation	Produced in rabbits immunized with E. coli-derived Human CALML5 fragment, and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human CALML5
Formulation:	PBS, pH7.0 with 0.03% Proclin300
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	IHC-P
----------------------	-------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:500-1:2000
--------------	---------------------


Please Note: Optimal concentrations/dilutions should be determined by the end user.

CALML5 Antibody, Rabbit PAb, Antigen Affinity Purified




Sino Biological
Biological Solution Specialist

Catalog Number: 207968-T08



Immunohistochemical staining of human CALML5 in human skin with rabbit polyclonal antibody at 1:1000 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of human CALML5 in human breast with rabbit polyclonal antibody at 1:1000 dilution, formalin-fixed paraffin embedded sections.