

# CD79B Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological  
Biological Solution Specialist

Catalog Number: 29750-T16

## GENERAL INFORMATION

|              |                                                                                                                                                                                                                                                    |
|--------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Immunogen:   | Recombinant Human CD79B Protein (Catalog#29750-H08H)                                                                                                                                                                                               |
| Preparation  | Produced in rabbits immunized with purified, recombinant Human CD79B (rh CD79B; Catalog#29750-H08H; NP_000617.1; Met1-Asp159). CD79B specific IgG was purified by Human CD79B affinity chromatography.                                             |
| Ig Type:     | Rabbit IgG                                                                                                                                                                                                                                         |
| Specificity: | Human CD79B                                                                                                                                                                                                                                        |
| Formulation: | 0.2 µm filtered solution in PBS                                                                                                                                                                                                                    |
| Storage:     | This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles. |

## APPLICATIONS

|               |       |
|---------------|-------|
| Applications: | ELISA |
|---------------|-------|

## RECOMMENDED CONCENTRATION

|       |                       |
|-------|-----------------------|
| ELISA | ELISA: 1:5000-1:10000 |
|-------|-----------------------|

***Please Note: Optimal concentrations/dilutions should be determined by the end user.***