

FHIT Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 310358-T40

GENERAL INFORMATION

Immunogen:	E. coli-derived Human FHIT fragment
Preparation	Produced in rabbits immunized with E. coli-derived Human FHIT fragment, and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human FHIT
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	AP3Aase,FRA3B

APPLICATIONS

Applications:	WB, IHC-P
---------------	-----------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:100-1:500
Western Blot	WB: 1:500-1:1000

Please Note: Optimal concentrations/dilutions should be determined by the end user.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)


Website: <http://www.sinobiological.com>

FHIT Antibody, Rabbit PAb, Antigen Affinity Purified




Sino Biological
Biological Solution Specialist

Catalog Number: 310358-T40



Immunohistochemical staining of human FHIT in human rectum with rabbit polyclonal antibody at 1:200 dilution, formalin-fixed paraffin embedded sections.



Anti-FHIT rabbit polyclonal antibody at 1:500 dilution

Lane A: HL-60 Whole Cell Lysate

Lysates/proteins at 30 µg per lane.

Secondary

Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution.

Developed using the ECL technique.

Performed under reducing conditions.

Predicted band size: 17 kDa

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>