

GOLPH2 / GOLM1 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 310897-T36

GENERAL INFORMATION

Immunogen:	E. coli-derived Human GOLPH2 / GOLM1 fragment
Preparation	Produced in rabbits immunized with E. coli-derived Human GOLPH2 / GOLM1 fragment, and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human GOLPH2 / GOLM1
Formulation:	PBS, pH7.0 with 0.03% Proclin300
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.
Alternative Names:	GOLM1

APPLICATIONS

Applications:	WB, IP
---------------	--------

RECOMMENDED CONCENTRATION

Western Blot	WB: 1:500-1:2000
Immunoprecipitation	IP:1-5μL/mg of lysate

Please Note: Optimal concentrations/dilutions should be determined by the end user.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)


Website: <http://www.sinobiological.com>

GOLPH2 / GOLM1 Antibody, Rabbit PAb, Antigen Affinity Purified




Sino Biological
Biological Solution Specialist

Catalog Number: 310897-T36



Anti-GOLPH2/GOLM1 rabbit polyclonal antibody at 1:500 dilution Lane A: MDA-MB-231 Whole Cell Lysate Lysates/proteins at 30 µg per lane. Secondary Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution. Developed using the ECL technique. Performed under reducing conditions. Predicted band size: 45 kDa



GOLM1 was immunoprecipitated using: Lane A: 0.5 mg MDA-MB-231 Whole Cell Lysate 4 µL anti-GOLM1 rabbit polyclonal antibody and 60 µg of Immunomagnetic beads Protein A/G. Primary antibody: Anti-GOLM1 rabbit polyclonal antibody, at 1:100 dilution Secondary antibody: Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution Developed using the ECL technique. Performed under reducing conditions. Predicted band size: 45 kDa Observed band size :75 kDa

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>