

flagellin / flaA Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 40067-RP01

GENERAL INFORMATION	
Immunogen:	40067-V07E
Preparation	Produced in rabbits immunized with purified, recombinant <i>Listeria monocytogenes</i> flagellin / flaA (r <i>Listeria monocytogenes</i> flagellin / flaA; Catalog#40067-V07E; YP_013330.1; Met1-Ser287). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	<i>Listeria monocytogenes</i> flagellin / flaA
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.
Alternative Names:	FlaA
APPLICATIONS	
Applications:	ELISA
	(Antibody's applications have not been validated with corresponding viruses. Optimal concentrations/dilutions should be determined by the end user.)
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:1000-1:2000

Please Note: Optimal concentrations/dilutions should be determined by the end user.